



*Sepax Technologies, Inc.*

5 Innovation Way, Newark, DE 19711, USA  
Tel: (302)366-1101. Fax: (302)366-1151.  
Toll-free: 1-877-SEPAX-US  
www.sepax-tech.com

---

## **Sepax Monomix HC/MC Ion Exchange Resins**

# **User Manual**

### Content

- I. Product Introduction
- II. Instructions for Use
- III. Storage
- IV. Ordering Information




## I. Product Introduction

Sepax Monomix HC/MC chromatographic resins are specially designed for analysis and purification of biological molecules. The resins are composed of hydrophilic polymethacrylate base beads with high physical strength and chemical stability. The monosized resin has an average particle size of 15, 30 and 60  $\mu\text{m}$  and an average pore size of 500 and 1000  $\text{\AA}$ . The resin surface is highly hydrophilic, which minimizes non-specific binding with biological samples. On the hydrophilic resin surface, strong cation exchange group (SP), strong anion exchange group (Q) and weak anion exchange group (DEAE) are covalently attached via proprietary technologies developed by Sepax. High mechanic stability ensures the resins tolerate high-pressure operation up to 10 bars, allowing for faster flow rate and shorter operation time. High chemical stability ensures the resins stand the acoustic cleaning procedure with a pH up to 14. These superior characters make the Sepax Monomix HC/MC IEX resins uniquely applicable at various bio-separation stages, from laboratory discovery, pilot-scale purification to industrial process.

Sepax Monomix HC/MC IEX resins are made with different linkers between the resin surface and the IEX functional groups. HC series are of much higher dynamic capacity than the MC series. The unique structural difference between these two types of resins offer different selectivity in the separation of biomolecules, such as peptides, insulins, proteins, antibodies, viruses, vaccines, heparins, and nucleotides.

### Resin Structure



**Figure 1.** Monomix HC/MC IEX resin structures. Strong cation exchange functional group is sulfonic acid (SP). Strong anion exchange functional group is trimethyl quaternary amine (Q), and weak anion exchange group functional group is diethyl amine (DEAE).



### Resin Characteristics

- High binding capacity and excellent biocompatibility
- Capable of withstanding high pressure and high flow rate due to rigid base bead matrix
- High separation resolution, efficiency and recovery
- High lot-to-lot consistency
- Easy to scale up
- Highly hydrophilic surface with minimal non-specific binding
- Small volume change under normal packing conditions
- Supply capacity: >100L per order

### Technical Specifications

**Table 1.** Monomix HC IEX Resins

| Resin Type                             | Monomix HC-Q   |            |            | Monomix HC-DEAE                                 |            |            | Monomix HC-SP      |                 |                 |
|--|--|------------|------------|---|------------|------------|--------------------|-----------------|-----------------|
| <b>Matrix</b>                          | Hydrophilic polymethacrylate   |            |            |   |            |            |                    |                 |                 |
| <b>Functional Group</b>                | -N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub>  |            |            | -N(C <sub>2</sub> H <sub>5</sub> ) <sub>2</sub> |            |            | -SO <sub>3</sub> H |                 |                 |
| <b>Particle Size (µm)</b>              | 15   | 30         | 60         | 15  | 30         | 60         | 15                 | 30              | 60              |
| <b>Pore Size (Å)</b>                   | 1000   |            |            | 1000  |            |            | 1000               |                 |                 |
| <b>DBC* (per mL resin)</b>             | ≥ 95mg BSA   | ≥ 90mg BSA | ≥ 80mg BSA | ≥ 95mg BSA                                      | ≥ 90mg BSA | ≥ 80mg BSA | ≥ 110mg Lysozyme   | ≥ 90mg Lysozyme | ≥ 90mg Lysozyme |
| <b>Maximum Linear Flow rate (cm/h)</b> | 1000   |            |            |   |            |            |                    |                 |                 |
| <b>Operation Temperature (°C)</b>      | ≤ 40   |            |            |   |            |            |                    |                 |                 |
| <b>pH range</b>                        | 2-12   |            |            |   |            |            |                    |                 |                 |
| <b>Maximum Pressure</b>                | ≤ 1 MPa (10 bar)   |            |            |   |            |            |                    |                 |                 |
| <b>Compatible Mobile Phases</b>        | Compatible with aqueous solution, a mixture of water and acetonitrile, ethanol, etc.<br>Typical buffers: Tris, phosphate, and acetate. |            |            |   |            |            |                    |                 |                 |
| <b>Storage</b>                         | 70% (v/v) in 20% ethanol   |            |            |   |            |            |                    |                 |                 |
| <b>Regeneration</b>                    | 1-2 M NaCl   |            |            |   |            |            |                    |                 |                 |
| <b>CIP</b>                             | 0.5 M HCl or 0.5-1.0 M NaOH. Impurities with strong hydrophobic binding can be cleaned with 0.1-1% Tween or Triton X-100               |            |            |   |            |            |                    |                 |                 |



**Table 2. Monomix MC IEX Resins**

| Resin Type                             | Monomix MC-Q   |            |            | Monomix MC-DEAE                                 |            |            | Monomix MC-SP      |                 |                 |
|--|--|------------|------------|---|------------|------------|--------------------|-----------------|-----------------|
| <b>Matrix</b>                          | Hydrophilic polymethacrylate   |            |            |   |            |            |                    |                 |                 |
| <b>Functional Group</b>                | -N <sup>+</sup> (CH <sub>3</sub> ) <sub>3</sub>  |            |            | -N(C <sub>2</sub> H <sub>5</sub> ) <sub>2</sub> |            |            | -SO <sub>3</sub> H |                 |                 |
| <b>Particle Size (µm)</b>              | 15   | 30         | 60         | 15  | 30         | 60         | 15                 | 30              | 60              |
| <b>Pore Size (Å)</b>                   | 500  |            |            | 500   |            |            | 500                |                 |                 |
| <b>DBC* (per mL resin)</b>             | ≥ 50mg BSA   | ≥ 40mg BSA | ≥ 30mg BSA | ≥ 50mg BSA                                      | ≥ 40mg BSA | ≥ 30mg BSA | ≥ 60mg Lysozyme    | ≥ 40mg Lysozyme | ≥ 30mg Lysozyme |
| <b>Maximum Linear Flow rate (cm/h)</b> | 1000   |            |            |   |            |            |                    |                 |                 |
| <b>Operation Temperature (°C)</b>      | ≤ 40   |            |            |   |            |            |                    |                 |                 |
| <b>pH range</b>                        | 2-12   |            |            |   |            |            |                    |                 |                 |
| <b>Maximum Pressure</b>                | ≤ 1 MPa (10 bar)   |            |            |   |            |            |                    |                 |                 |
| <b>Compatible Mobile Phases</b>        | Compatible with aqueous solution, a mixture of water and acetonitrile, ethanol, etc.<br>Typical buffers: Tris, phosphate, and acetate. |            |            |   |            |            |                    |                 |                 |
| <b>Storage</b>                         | 70% (v/v) in 20% ethanol   |            |            |   |            |            |                    |                 |                 |
| <b>Regeneration</b>                    | 1-2 M NaCl   |            |            |   |            |            |                    |                 |                 |
| <b>CIP</b>                             | 0.5 M HCl or 0.5-1.0 M NaOH. Impurities with strong hydrophobic binding can be cleaned with 0.1-1% Tween or Triton X-100               |            |            |   |            |            |                    |                 |                 |

\* Dynamic Binding Capacity (DBC) measurement conditions:

- HC/MC-Q: linear flow rate 180 cm/h, 2.0 mg/mL of BSA in 50 mM Tris buffer, pH 8.5
- HC/MC-DEAE: linear flow rate 180 cm/h, 2.0 mg/mL of BSA in 20 mM Tris buffer, pH 8.0
- HC/MC-SP: linear flow rate 360 cm/h, 1.0 mg/mL of lysozyme in 50 mM phosphate buffer, pH 6.0

## II. Instructions for Use

### 2.1 Safety Precautions

For information on the safe use of this product, please refer to the Safety Data Sheet (SDS).

### 2.2 Resin Preparation before Use

This product is normally shipped in an aqueous solution containing 20% ethanol and should be washed with deionized water before getting in contact with salt buffer



solutions. This can be accomplished by rinsing with 3 times column volume of deionized water and can be part of the column packing procedure (see Section 2.3.2 below).

### 2.3 Column Packing

- 2.3.1 Calculate column volume (CV):  $CV = \text{column cross-sectional area } (\pi r^2) \times \text{bed height } (h)$ ,  $r$  is the column radius;
- 2.3.2 Gently stir the resin to completely disperse it to form a uniform slurry. Measure the required volume of the resin slurry and pour it into a clean transparent glass or plastic vessel. After natural sedimentation, decant 20% ethanol solution. Add 3 CV of deionized water, gently stir evenly and then settle for about 30 minutes, decant the supernatant, and repeat 3 times;
- 2.3.3 After removing the supernatant, add column packing buffer (0.5M NaCl solution) to make 60-70% resin slurry (volume based), stir evenly and soak the resin for more than 12 hours (overnight);
- 2.3.4 Measure about 1.1 times CV of the resin and gently stir evenly to make the final packing slurry. Pour it into a column with filter plate at the bottom to allow packing buffer to flow through and resin to settle steadily;
- 2.3.5 Place a flow distributor on the top of the column, press down the resin bed and connect to a pump;
- 2.3.6 Use 2-3 CV column packing buffer to flush the resin bed, compress resin beds at 2 times of normal working flow rate. The distributor position can be adjusted during the compressing process to ensure tightness of the resin bed. It is not recommended to use suction or gravity-only sedimentation to pack a column, especially for columns with a bed height of more than 10 cm;
- 2.3.7 Evaluation of column packing quality is carried out using a low molecular weight, unretained compound. The specific operating parameters are as follows:

| Sample         | 1.0% (v/v) acetone aqueous solution   | 0.8 M NaCl                   |
|----------------|---|------------------------------|
| Sample Volume  | 1.0–2.0% of bed volume  | 1.0–2.0% of bed volume       |
| Mobile Phase   | Water or diluted buffer   | Water or 0.5 M NaCl solution |
| Flow Rate      | 180 cm/h  | 180 cm/h                     |
| Detector       | UV 280 nm   | Conductivity                 |
| Specifications | Tailing factor: 0.8 – 1.5<br>30µm resin column efficiency: $\geq 4,000/m$<br>60µm resin column efficiency: $\geq 2,000/m$ |                              |

- 2.3.8 In case of non-ideal result, such as peak tailing, solutions include:
  - Reduce the concentration of slurry



---

- Increase packing flow rate

- Extend packing time

If peak fronting occurs, solutions are the opposite of the above.

## **2.4 Column Use**

- 2.4.1 Select and optimize a buffer system according to specific characteristics (isoelectric point, ion-exchange group, etc.) of targeted molecules to be separated/purified or analyzed;
- 2.4.2 Equilibrate the column with about 5 CV of buffer until conductivity and pH of effluent are constant, consistent with those of the fresh buffer;
- 2.4.3 **Sample preparation:** Solid samples can be dissolved in a buffer solution; low concentration samples can be dialyzed with a buffer to increase concentration; high concentration samples can be diluted with a buffer. Samples with insoluble impurities should be filtered first to avoid clogging the column and to prolong column life;
- 2.4.4 **Loading:** A sample loading volume should be determined according to resin loading capacity and purity of targeted molecule in crude sample; after loading is completed, continue to pump buffer until a stable baseline is obtained;
- 2.4.5 **Elution:** According to the characteristics of targeted molecule, choose a method (increasing the salt concentration or changing pH) to elute targeted molecule bound to the resin bed;
- 2.4.6 **Regeneration:** After each use of purification/separation, rinse the column with 1-2M NaCl solution to remove impurities adsorbed on the resin bed;
- 2.4.7 **Cleaning-in-Place (CIP):** If impurities cannot be removed through the regeneration step, causing column clogging, increased back pressure, or decreased flow rate, the performance of the column can be restored by forward or reverse online cleaning – normally washed with 5 CV of 0.5 M NaOH. In general, online cleaning can result in increased back pressure of the column, so it is recommended to use linear flow rates within 0.5 times of normal operation conditions. Specific CIP methods are determined based on the characteristics of the impurities, some examples below:
- For impurities such as precipitated or denatured substances, wash with 5 CV of 0.5-1.0 M NaOH; if the result is not satisfactory, use 0.5-1.0 M NaOH plus 1M NaCl, and then use at least 5 CV of 0.22um filtered buffer (pH6-8) to wash the column;
  - For impurities bound with strong hydrophobicity, wash the column with 2 CV of non-ionic detergent (such as Tween or Triton X-100 at a concentration of 0.1-1%), and then immediately use at least 5 CV of filtered buffer (pH6-8) to wash the column; it can also be washed with 3-4



## Sepax Technologies, Inc.

5 Innovation Way, Newark, DE 19711, USA  
Tel: (302)366-1101. Fax: (302)366-1151.  
Toll-free: 1-877-SEPAX-US  
www.sepax-tech.com

CV of 70% ethanol or 30% isopropanol, and then at least 5 CV of filtered buffer (pH6-8) to wash.

### III. Storage

Chromatographic resins that will not be immediately used should be stored in an aqueous solution containing 20% ethanol at 4 to 35 °C in a sealed container; resins that have packed in a column can be stored in an aqueous solution containing 20% ethanol at 4 to 35 °C.

### IV. Ordering Information

| Resin             | Particle Size, Pore Size | Product Number | Package Size (L)   | Prepacked Column (mL) |
|-------------------|--------------------------|----------------|--------------------|-----------------------|
| Monomix HC15-Q    | 15µm, 1000Å              | 280815950      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix HC15-DEAE | 15µm, 1000Å              | 280915950      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix HC15-SP   | 15µm, 1000Å              | 280615950      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix MC15-Q    | 15µm, 500Å               | 280415500      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix MC15-DEAE | 15µm, 500Å               | 280515500      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix MC15-SP   | 15µm, 500Å               | 280215500      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix HC30-Q    | 30µm, 1000Å              | 280830950      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix HC30-DEAE | 30µm, 1000Å              | 280930950      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix HC30-SP   | 30µm, 1000Å              | 280630950      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix MC30-Q    | 30µm, 500Å               | 280430500      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix MC30-DEAE | 30µm, 500Å               | 280530500      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix MC30-SP   | 30µm, 500Å               | 280230500      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix HC60-Q    | 60µm, 1000Å              | 280860950      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix HC60-DEAE | 60µm, 1000Å              | 280960950      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix HC60-SP   | 60µm, 1000Å              | 280660950      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix MC60-Q    | 60µm, 500Å               | 280460500      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix MC60-DEAE | 60µm, 500Å               | 280560500      | 0.5, 1, 5, 10, 100 | 1, 5                  |
| Monomix MC60-SP   | 60µm, 500Å               | 280260500      | 0.5, 1, 5, 10, 100 | 1, 5                  |

For resin products or prepacked columns not listed above, please contact us.