

Chromatography Column Recommendation for Corresponding Sepax Brand

	Sepax	Tosoh	Waters	Agilent	Thermo	Phenomenex
SEC	SRT SEC-100 SRT-C SEC-100	—	—	Agilent Bio SEC-5 100Å	BioBasic SEC 120	—
	SRT SEC-150 SRT-C SEC-150	TSKgel G2000SWxL	—	Agilent Bio SEC-5 150Å	—	BioSep-SEC-S2000
	SRT SEC-300 SRT-C SEC-300	TSKgel G3000SWxL	—	Agilent Bio SEC-5 300Å	BioBasic SEC 300 MAbPac SEC-1	BioSep-SEC-S3000
	SRT SEC-500 SRT-C SEC-500	TSKgel G4000SWxL	Xbridge Protein BEH 450 SEC	Agilent Bio SEC-5 500Å	—	BioSep-SEC-S4000
	SRT SEC-1000 SRT-C SEC-1000	—	—	Agilent Bio SEC-5 1000Å	BioBasic SEC 1000	—
	SRT SEC-2000 SRT-C SEC-2000	—	—	Agilent Bio SEC-5 2000Å	—	—
	Zenix SEC-80	—	—	—	—	—
	Zenix SEC-100	TSKgel Super SW G2000	Xbridge Protein BEH 125 SEC	Agilent Bio SEC-3 100Å	—	—
	Zenix SEC-150	TSKgel Super SW G2000	Xbridge Protein BEH 125 SEC	Agilent Bio SEC-3 150Å	—	—
	Zenix SEC-300	TSKgel Super SW G3000	Xbridge Protein BEH 200 SEC BioResolve SEC mAb Column	Agilent Bio SEC-3 300Å Agilent AdvanceBio SEC300Å	—	Phenomenex Yarra SEC 3000
Unix-C SEC-300	TSKgel UP-SW3000	Acquity UPLC Protein BEH 200 SEC	—	—	bioZen 1.8um SEC-3	
IEX	Proteomix SCX	TSKgel SP-5PW	Protein-Pak Hi Res IEX SP BioResolve SCX mAb	Agilent Bio IEX SCX	ProPac SCX-10 MAbPac SCX	—
	Proteomix WCX	TSKgel CM-5PW	Protein-Pak Hi Res IEX CM	Agilent Bio IEX WCX	ProPac WCX-10 ProPac Elite WCX	bioZen WCX
	Proteomix SAX	TSKgel SuperQ-5PW TSKgel QAE-2SW	Protein-Pak Hi Res IEX Q	Agilent Bio IEX SAX	ProPac SAX-10	—
	Proteomix WAX	TSKgel DEAE-5PW	—	Agilent Bio IEX WAX	ProPac WAX-10	—
	Antibodix	TSKgel CM-STAT	Protein-Pak Hi Res CM	Agilent Bio IEX WCX	ProPac WCX-10	bioZen WCX
HIC	Proteomix HIC Phenyl	TSKgel Phenyl-5PW	—	—	—	—
	Proteomix HIC Butyl	TSKgel Butyl-NPR	Protein-Pak Hi Res HIC	AdvanceBio HIC	ProPac HIC-10 MAbPac HIC-10	—
	Proteomix HIC Propyl	—	—	—	—	—
	Proteomix HIC Ethyl	TSKgel Ether-5PW	—	—	—	—
AFC	ProAqa Excel	TSKgel ProteinA-5PW	—	Agilent Bio-Monolith Protein A	MAbPac Protein A POROS proteinA	—



Chromatography Column Recommendation for Corresponding Sepax Brand

	Sepax	Bio-Rad	Waters	Thermo	Phenomenex
Sugar	Carbomix H-NP	Aminex HPX-87H Fast Acid	Sugar-Pak 1	HyperRez XP Carbohydrate H+	Rezex RHM- Monosaccharide RFQ-Fast Acid
	Carbomix Ca-NP	HPX-87Ca	—	HyperRez XP Carbohydrate Ca+	RCU-USP Sugar Alcohols
	Carbomix Pb-NP	HPX-87P	—	HyperRez XP Carbohydrate Pb+	RezexRHM- Monosaccharide
	Carbomix K-NP	HPX-87K	—	—	Rezex RKP- Potassium
	Carbomix Na-NP	HPX-87N	—	HyperRez XP Carbohydrate Na+	RNO- Oligosaccharide

ADDITIONAL INFORMATION

A creation of a customer account will provide you with:

- Full access to application notes of your interested sample type
- Method development training webinars
- Up-to-date listing of where our columns have been cited in the scientific literature
- Easy way to view prices and request quotes



 www.sepax-tech.com

 info@sepax-tech.com