

Sepax Technologies, Inc. - Column Accessories

ANALYTICAL
Guard

Filter

ANALYTICAL

For Analytical Columns
with Particle Size 5 µm and 10 µm

- PN: [102000-P355](#) PEEK Precolumn Filter & Frit (2 µm)
- PN: [102001-P355](#) PEEK Refill Frits (2 µm), 5 units/pk
- PN: [102000-S355](#) SS Precolumn Filter & Frit (2 µm)
- PN: [102001-S355](#) SS Refill Frits (2 µm), 5 units/pk

with Particle Size below 3 µm


- PN: [102000-P356](#) PEEK Precolumn Filter & Frit (0.5 µm)
- PN: [102001-P356](#) PEEK Refill Frits (0.5 µm), 5 units/pk
- PN: [102000-P346](#) PEEK Precolumn Filter SS Frit (0.5 µm)
- PN: [102001-P346](#) SS Refill Frits (0.5 µm), 5 units/pk



PEEK
Column
Coupler

[102002-COUPLER](#)

Analytical
PEEK
Coupler
with
Flexible
Tubing



[102006-COUPLER](#)

Cartridge & Holder

Packed with Resin,
Part Number is dependent to Phase
[Contact sales for assistance.](#)

Replacement Cartridge Holder

- PN: [102000-2001](#)
Holder for 2.1mm ID column
- PN: [102000-4001](#)
Holder for 4.6 mm ID column



PREP
Filter

For Prep Column (≥ 21.2mm ID)
Particle Size 5 µm or Above

- PN: [102020-21200](#)
Stainless Steel Precolumn Filter & Frits (2 µm)
- PN: [102020-00001](#)
Stainless Steel Refill Frits, (2 µm),
5 units/pk



PEEK/SS
Coupler



[102003-COUPLER](#)

Cartridge & Holder

Packed with Resin,
Part Number is dependent to Phase
[Contact Sales for assistance.](#)



PREP
Guard